



# BD® Research Cloud with BD Horizon™ Panel Maker

Accelerate your workflow with BD® Research Cloud (BD RC) panel creation module, now featuring BD Horizon™ Panel Maker, an intelligent solution for automated panel design. Enjoy seamless integration with instrument configurations, strategic reagent recommendations and scientifically validated resources; all crafted to save time and deliver high-performing panels.



Note: All panels should be validated experimentally based on your specific experimental conditions.

## Instrument selection

The screenshot shows the BD FACSDiscover A8 (Spectral) software interface. At the top, it displays 'Cytometer Summary' for 'BD FACSDiscover A8 (Spectral)' with 'CS1\* Board Lot File'. Below this, 'CONFIGURATION NAME' is set to 'SL (22)IV\_3DV\_148.12Y\_0\_8R\_1Tmg (38)'. The main panel shows 'Power (mW): 60' and 'Wavelength: 349'. The 'Filters' section lists various filter combinations with their corresponding 'Parameter Name(s)'. Below this is a list of four BD instruments: BD FACSDISCOVER A8, BD FACSDISCOVER S8, BD FACSYMPHONY S6 SE, and BD FACSYMPHONY A5 SE, each with a small image and a green circular icon.

## Reagent optimization

The screenshot shows the BD Horizon Panel Maker software interface. It features a 'Spectral Hotspot Matrix' on the left with a color-coded legend for reagent concentration (Low, Mid, High, Unknown, Critical) and a 'Panel Maker Summary' table on the right. The table lists various markers and their panel assignments across different instruments.

## Antigen population guide

The screenshot shows the BD Horizon Panel Maker software interface. It displays an 'Antigen Population Guide' table. The columns are 'NAME', 'Tregs', 'CD4 naïve T-cells', 'CD4 memory T-cells', 'CD8 naïve T-cells', and 'CD8 memory T-cells'. The rows list various antigens: CD3, CD4, CD8, CD45RA, CD127 (IL-7 Receptor α), CD25 (IL-2 Receptor α), and CD197 (CCR7). Each cell in the table contains a colored dot representing the phenotype for that specific antigen and cell type combination.

## Instrument optical performance

- Get optimized panel recommendations driven by preloaded, measured data for enhanced accuracy. Cytometer-specific\* spectral signature data generated using fluorochrome-stained cells.

\*BD Horizon™ Panel Maker is available for BD FACSDiscover™ S8 Cell Sorter, BD FACSDiscover™ A8 Cell Analyzer, BD FACSymphony™ A5 SE Cell Analyzer, BD FACSymphony™ S6 SE Cell Sorter and Cytek Aurora flow cytometers. Measured data was generated using surface-stained CD4 and human PBMCs.

- Achieve spectral results that closely reflect real-world performance by leveraging your instrument's optical configuration.
- Simplify complex workflows with the only panel designer that supports integrated imaging and spectral panel design for BD FACSDiscover™ Instruments.


## Strategic reagent recommendations

- Get up to three optimized panel suggestions quickly based on your selected markers and inputs.
- Rapidly compare and evaluate panels using visualization tools like the BD® Spectral Hotspot Matrix.
- Streamline your workflow with direct access to BD's extensive reagent catalog.
- Enjoy flexibility and choice with the ability to select dyes from multiple manufacturers.

## Scientific resources

- Design with confidence whether starting from scratch or using curated databases of published panels and populations phenotypes.
- Leverage BD Horizon™ Panel Maker design rules and tools backed by decades of BD expertise for reliable results.
- Leverage high-quality scientific insights throughout the panel design process.

# Guided workflow



## 1 Select a cytometer

Cytometer-specific, preloaded measured data is combined with instrument optical configurations for results that closely mirror real-world performance.

## 2 Select targets and define populations

Input key biology and experiment details: imaging marker selection, clones, antigen density, and exclusion markers.

## 3 Create Panel

### BD Horizon™ Panel Maker




(Automated Panel Design)

Use the BD Horizon™ Panel Maker automation to pair your biological and experimental inputs with ideal fluorochrome combinations. Then use powerful imaging and spectral tools like BD Spectral Hotspot Matrix, heat maps and spectral graphs to visualize panel recommendations and evaluate their performance.

## 4 Select reagents and fluorochromes

View and choose reagents from your current inventory, BD's integrated catalog, your custom catalog and dyes from multiple manufacturers.



## 5 Create shopping list

Generate quotes and purchase directly through BD's eCommerce platform or send a shopping list to your procurement team.

Open Cart at [bdbiosciences.com](http://bdbiosciences.com) for special pricing


[Add to BD Cart](#)

Open Shopping List at [bdbiosciences.com](http://bdbiosciences.com)

[Add to Shopping List](#)

# Leverage proven best practices with BD Horizon™ Panel Maker optimization

Capitalizing on over 50 years of BD's panel design expertise, the BD Horizon™ Panel Maker algorithm incorporates our experts' tools, rules and best practices into its framework.



## Resources

### FlowJo University

[flowjo.com/learn/flowjo-university/bd-research-cloud](http://flowjo.com/learn/flowjo-university/bd-research-cloud)

### Webinars

[flowjo.com/learn/webinars](http://flowjo.com/learn/webinars)

### Documents

[flowjo.com/learn/marketing-materials](http://flowjo.com/learn/marketing-materials)

**Sign-up for free:** [bdresearchcloud.com/login](http://bdresearchcloud.com/login)

For Research Use Only. Not for use in diagnostic or therapeutic procedures.  
BD Life Sciences, Milpitas, CA 95035, U.S.

**bdbiosciences.com**

